

# DISINFECTANT SPRAY

Textiles & surfaces

**Hydroalcoholic disinfectant spray**, delicately perfumed, intended for quick disinfection of **textiles** and **any other surfaces**.

✓ **Active against enveloped viruses as Coronavirus, H1N1 (influenza A) and rotavirus.**

in compliance with EN 14476

✓ **Delicate fragrance**

✓ **Suitable for disinfection of textiles**

✓ **1 aerosol = 10m<sup>2</sup>**

Recommended within the **medical field** (hospitals, medical institutions or veterinary field), **communities, offices, public places** or **domestic**.

Suitable for disinfection of **textiles as seats, sofas, chairs, stretchers** etc.

Use adapted for disinfection in the field of cars and motorcycles : **inside the helmet, shoes** and **interiors of vehicles** (automotive seat, dashboard etc.).

Use on objects at risk of manual transmission : **door handles, phones, switches, computer equipment** (mouse, keyboard), etc.

⚠ For delicate surfaces, spray beforehand on a microfibre to a proper disinfection.

Disinfectant properties in compliance with standards at 20°C, clean condition :

NORMES	STRAINS	CONDITIONS
Bactericidal in compliance with EN 13727	Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Enterococcus hirae	In 30 secondes
Bactericidal in compliance with EN 13697		In 5 minutes
Fungicidal in compliance with EN 13624	Candida Albicans and Aspergillus brasiliensis	In 30 secondes
Fungicidal in compliance with EN 13697		In 15 minutes
Limited spectrum virucidal in compliance with EN14476	Murine norovirus (in 30 secondes), Adenovirus (in 30 secondes) et Vacciniavirus (in 2 minutes)	In 30 secondes and 2 minutes
	H1N1 (grippe A) and Rotavirus	In 30 secondes
Mycobactericidal in compliance with EN 14348	Mycobacterium avium complex and mycobacterium terrae (= tuberculocidal activity)	En 5 minutes



300 ML - Ref : A22472

[Technical sheet](#)

Click here to download